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The main content area is divided into two sections. The left section, titled '検査日時' (Examination Date), shows a list of examinations: 2021/10/01 11:06:56 (CT) and 2015/10/01 17:16:34 (CT, 胸部, 読影 檀山). The right section, titled '第1読影' (First Reading), shows the reading process. It includes a 'RF判定' (RF Judgment) section with a 'なし' (None) result. Below this, there are three reading slots (第1所見, 第2所見, 第3所見) with dropdown menus for '左下葉' (Left Lower Lobe), '陳旧性結核' (Old Tuberculosis), and 'C. 一年後に低線' (C. Low line in one year). To the right of these slots are three CT scan images labeled '第1所見', '第2所見', and '第3所見'. Below the reading slots is a 'コメント' (Comment) field. The '第2読影' (Second Reading) section follows a similar layout, with a 'mens onboarding1' reading doctor and a '左肺下葉S9に小結節影が見られますが、経時的変化はありません。' (A small nodule shadow is seen in the S9 segment of the left lower lung, but there is no change over time.) comment. At the bottom, there's a '総合所見' (Overall Finding) section with a text area containing '軽度の異常が見られますが治療は必要ありません。来年度も低線量CT検診を受診してください。' (Mild abnormalities are seen, but treatment is not necessary. Please undergo low-dose CT examination next year.) and a '保存' (Save) button. The footer contains navigation buttons like '定型文' (Template Text), '一次確定' (First Confirmation), '二次確定' (Second Confirmation), '印刷' (Print), '改訂' (Revision), 'CAD参照' (CAD Reference), and '変更履歴' (Change History).

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